



- MS238-A 4.4mm** 65 Dia. 44x.005, 6x.008, 4x.013, 6x.02, 2x.025, 2x.037, 1x.50= 1.06 ct. 6.0 Gr.  
(1, 1.2, 1.4, 1.6, 1.8, 2.1, 4.4mm cushion)
- MS238-A 5.5mm** 69 Dia. 48x.005, 6x.008, 4x.013, 6x.02, 2x.025, 2x.037, 1x1.00= 1.58 ct. 6.1 Gr.  
(1, 1.2, 1.4, 1.6, 1.8, 2.1, 5.5mm cushion)
- MS238-A 6.0mm** 71 Dia. 50x.005, 6x.008, 4x.013, 6x.02, 2x.025, 2x.037, 1x1.25= 1.84 ct. 6.2 Gr.  
(1, 1.2, 1.4, 1.6, 1.8, 2.1, 6.0mm cushion)
- MS238-A 6.5mm** 75 Dia. 54x.005, 6x.008, 4x.013, 6x.02, 2x.025, 2x.037, 1x1.50= 2.11 ct. 6.4 Gr.  
(1, 1.2, 1.4, 1.6, 1.8, 2.1, 6.5mm cushion)
- MS238-A 7.0mm** 77 Dia. 56x.005, 6x.008, 4x.013, 6x.02, 2x.025, 2x.037, 1x2.00= 2.62 ct. 6.5 Gr.  
(1, 1.2, 1.4, 1.6, 1.8, 2.1, 7.0mm cushion)